

No foot control needed.



SCALEX® Glide

HAND-ACTIVATED PIEZO ULTRASONIC SCALER

- Revolutionary hand-activated piezo ultrasonic scaler, U.S. and world patents pending.
- No foot control needed.
- Futuristic touch sensitive panel technology.
- Autoclavable handpiece with touch ring is accessible from 360°.
- State-of-the-art piezoelectric technology is self-tuning and load-adjusting for smooth and effortless operation.
- Advanced water-cooled circuitry designed for greater stability during operation.
- Extremely gentle on the patient.
- Wide selection of high quality tips to correspond with specific clinical applications.
- Available for 100-240V 50/60Hz operation.

IT'S ALL IN THE HANDS.

DENTAMERICA®

VALUE BY DESIGN

360° Trigger Ring



Autoclavable Handpiece

PIEZO SCALER TIPS



88002

Reorder No.155-03
Endo Chuck



88003

Reorder No.155-01
Universal A



88004

Reorder No.155-51
Universal B



88005

Reorder No.155-02
Chisel



88006

Reorder No.155-48
Perio

SPECIFICATIONS

INPUT: 100V~240V 50/60Hz

OUTPUT: DC 24V 1.5A

FREQUENCY: 27,000 to 32,000 Hz

DIMENSIONS:
W221mm x D142mm x H43mm

WEIGHT: 620g

DIMENSIONS(HP): Ø22mmxW110mm

WEIGHT(HP): 50g

ISO 13485 **GMP**

DENTAMERICA INC.

18688 E. San Jose Ave. City of Industry, CA 91748 U.S.A.
TEL: (626) 912-1388 FAX: (626) 913-0510
www.dentamerica.com
www.dentamerica.asia

DISTRIBUTOR

DENTAMERICA®