

# ROTEX™ 782E

## COMPACT ELECTRIC E TYPE HANDPIECE UNIT



ISO 9001  
ISO 13485



GMP

- Ball bearing high torque electric motor.
- 1,000 to 35,000 R.P.M. variable speed.
- Double insulation heavy duty construction.
- Automatic cooling fan for continuous high speed performance.
- Forward/Reverse switch.
- Application for grinding, polishing, milling of precious metals, ceramics, acrylics and steel.
- Standard package includes Variable Speed Foot Control, Handpiece and Handpiece Rest. Circuit overload protection.
- Unit is available for 100/110V 50/60Hz or 220/240V 50Hz operation.

ROTEX™  
782E

**DENTAMERICA®**

# ROTEX™ 782E

## COMPACT ELECTRIC E TYPE HANDPIECE UNIT

Reorder No. 7822

Straight Attachment



Reorder No. 7824

Contra Attachment


Reorder No. 773E

E-Type Motor Handpiece



Reorder No. 775

Variable Speed Foot Control



### SPECIFICATIONS :

- INPUT : 115V/230V 50/60Hz
- OUTPUT : DC 2V-39V
- SPEED : 1,000-35,000 rpm
- DIMENSIONS : L120xW190xH80mm
- WEIGHT : 2.0kg

### DENTAMERICA®

18688 E. San Jose Ave.  
City of Industry, CA 91748 U.S.A.  
Phone: (626) 912-1388  
Fax: (626) 913-0510  
www.dentamerica.com

DISTRIBUTOR

DENTAMERICA®